


HCV

One step qualitative immunochromatographic assay for detection of antibodies against HCV in human serum, plasma and whole blood


- Double Antigen Sandwich site Immunoassay
- Coated Ag: Core NS3, NS4 & NS5
- Specimen Type: Serum, Plasma and Whole Blood
- Specimen Volume: 10 μ l Serum / Plasma or 20 μ l Whole Blood
- Buffer Volume: 3 drops
- Result Interpretation Time: 20 Minutes
- Diagnostic Sensitivity: > 99.5 %
- Diagnostic Specificity: > 99.0 %
- Long shelf life



Test protocol



Interpretation of test results



Product Code	Material Description	Pack Size
RPDHCV - 01	MERISCREEN HCV	50 Test



Meril Diagnostics Pvt. Ltd.

D1-D3, Second Floor, Meril Park, Muktanand Marg, Chala, Vapi - 396 191, Gujarat, India.

T: +91 260 3052 100 F: +91 260 3052 125 E: diagnostics@merillife.com W: www.merillife.com