

# HIV 1-2 WB


One step qualitative immunochromatographic assay for detection of antibodies against HIV - 1 including subtype C and HIV - 2 in human serum, plasma and whole blood

- 3<sup>rd</sup> Generation Assay
- Detection of IgM, IgG & IgA Antibodies
- Coated Ag: gp120, gp41 for HIV 1 and gp36 for HIV 2
- Specimen Type: Serum, Plasma, Whole blood
- Specimen Volume: 10  $\mu$ l Serum / Plasma or 20  $\mu$ l Whole Blood
- Buffer Volume: 3 drops
- Result Interpretation Time: 20 Minutes
- Diagnostic Sensitivity: > 99.5 %
- Diagnostic Specificity: > 99.0 %
- Long shelf life




Evaluated by NARI | Procured by NACO

## Test protocol



## Interpretation of test results



Product Code	Material Description	Pack Size
RPDHIV - 01	MERISCREEN HIV 1-2 WB	50 Test